










サルファ剤  
Sulfonamides 磺胺药物

O950302C



- 1   
Sulfamethazine 磺胺二甲嘧啶
- 2   
Sulfadoxine 磺胺多辛
- 3   
Sulfamonomethoxine hydrate  
磺胺间甲氧嘧啶水合物
- 4   
Sulfisoxazole 磺胺二甲异唑
- 5   
Sulfamethoxazole 磺胺甲噁唑
- 6   
Sulfachloropyridazine  
磺胺氯达嗪
- 7   
Sulfadimethoxine 磺胺地索辛
- 8   
Sulfaquinoxaline sodium  
磺胺喹啉钠

Column : YMC-Pack ODS-A (5  $\mu$ m, 12nm)  
250  $\times$  4.6 mm I.D.  
Eluent : THF/water/acetic acid (30/70/1)  
Flow rate : 0.5 mL/min  
Temperature : 37  $^{\circ}$ C  
Detection : UV at 275 nm, 0.16 AUFS  
Injection : 5  $\mu$ L (42  $\mu$ g/mL)